



Hillrom™

Likorall® Overhead Lifts
200, 242, 250

Likorall® Overhead Lift

A THOUGHTFUL, USER-FRIENDLY DESIGN WITH A STRONG FOCUS ON SAFETY.

The Likorall lift is known for all these and more, making this overhead lift system a good fit for a variety of healthcare settings.

Safety First, Every Time

- Mechanical or electrical emergency lowering available
- Patented safety drum adds extra security and ensures even, controlled lowering
- Safety Squeeze Protection Limit Switch protects the lift motor from mechanical strain and helps prevent injuries to patients and caregivers
- Room-to-room (R2R) model helps you efficiently move patients between rooms without a continuous rail connection
- On the ES model, the hand control can help you navigate through switches and turn tables.
- UltraTwin ceiling mounted lift system consists of two Likorall lift motors that provide a combined lifting capacity.
- Safe Working Load
 - + Standard models: 200-250 kg / 440-550 lbs, depending on model
 - + UltraTwin model: 400-500 kg / 880-1,100 lbs, depending on model



Moving the Lift Motor (excluding LR242 S and LR 242 ES)

The Quick-Release Hook™ Carriage makes it easy to move the system between rooms.



In-Rail Charging (Likorall 250 S IRC only)

The lift unit is constantly charged along the complete rail system, so it's always ready for use.

Available Models

Likorall® Overhead Lift Systems 200 | 242 S, S R2R, ES, ES R2R
250 ES, S IRC



	LR200	LR242 S	LR242 S R2R	LR242 ES	LR242 ES R2R	LR250 ES	LR250 S IRC
Maximum Weight Capacity	200 kg / 440 lb	200 kg / 440 lb	200 kg / 440 lb	200 kg / 440 lb	200 kg / 440 lb	250 kg / 550 lb	250 kg / 550 lb
Lift Interval	2050 mm	2050 mm	2050 mm	2050 mm	2050 mm	2050 mm	2050 mm
Lift Speed	50 mm/s	50 mm/s	50 mm/s	50 mm/s	50 mm/s	50 mm/s	50 mm/s
Battery Technology	Lead-Acid	Lead-Acid	Lead-Acid	Lead-Acid	Lead-Acid	Lead-Acid	Lead-Acid
Quick Release Hook (W/O adapter)	X		X		X	X	X
Electrical Emergency Lowering	X	X	X	X	X	X	X
Manual Emergency Lowering		X		X			
Electrical Emergency Lowering Reachable From Floor	X						
Manual Emergency Lowering Reachable From Floor		X		X			
Room-to-Room			X		X		
Fixed Hook (W/O adapter)		X		X			
In-Rail Charging							X
Transfer Motor Compatibility				X		X	
Weight of Motor (kg)	12.5	13	13.2	13	13.2	12.6	12.6
Dimensions (L x W x H)	340 mm x 250 mm x 165 mm / 13.39" x 9.84" x 6.5"						



Hillrom™

Autorizované zastoupení pro ČR:

DARTiN

DARTIN, spol. s r.o.
Suchdolská 688, 252 62 Horoměřice
Tel./fax: 241 470 361, 261 216 076
www.dartin.cz
e-mail: info@artin.cz

Autorizované zastúpenie pre SR:

DARTiN

DARTIN Slovensko, s.r.o.
Ľ. Podjavorinskej 25, 080 05 Prešov
Tel./fax: 051/7734403, 7723634
www.dartin.sk
e-mail: info@artin.sk

ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at hillrom.com.

For more information, please contact your local distributor or Hillrom sales representative.

hillrom.com

OVERHEAD LIFTS

Likorall® overhead lift is a general-purpose lift with the intended use in: health care, intensive care and rehabilitation. Likorall overhead lift is designed and developed for the complete range of Hillrom fixed installed and free-standing lift systems.

Class 1

Manufacturer: Liko AB, Nedre vägen 100, 975 92 Luleå, Sweden

These products are regulated health care products which, where required by applicable regulations, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documentation accompanying medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.