



## A simple, reliable raising device.

Sabina- 200 sit-to-stand lift is used to assist people who have difficulty raising from a sitting to a standing position, for example from a wheelchair, chair, bed or when using the toilet.

The lift mast has three settings to match the height of the individual patient. The carefully designed accessories ensure that raising and standing practice can be adapted to suit the needs of the individual patient. This improves the conditions for safe and comfortable raising.

### Important Features

- A wide range of accessories for many individual adaptations.
- Adjustable height setting for different-sized patients.
- The lower leg support and the foot rest can easily be removed, when needed.

**The best way of raising yourself.**

Sabina™ is the result of comprehensive study and careful analysis of geometry and body movement patterns. Sabina's lifting function and anchor points close to the patient's centre of gravity make for the most effective way of raising yourself.

**Accessories for adaptation to the individual.**

There are four types of sit-to-stand vests available; select the one which suits your patient best. If the patient requires a little extra help to start the lift you can use a Seatstrap also. Sabina™ SideSupport provides extra security for patients that have a tendency to lean to the side.

**Accessories.**

- 2027101 Liko™ Sabina™ SideSupport
- 2027011 Heel support Sabina™
- 2027006 Seatstrap Slingbar
- 3591115 Liko™ Sabina™ Seatstrap

*Raising and base-width adjustment are easily controlled by the caregiver via the hand control.*



Further information on the various sit-to-stand vests Liko™ SupportVest™, Liko™ Solo™ SupportVest™, Liko™ SafetyVest™ and Liko™ SafetyVest™ with extra back support is available at [www.liko.com](http://www.liko.com)

- Max. load: 160 kg
- Emergency Lowering: Electrical
- Material: Powder-coated Steel
- Prod. number: 2020020



Sabina™ 200 has been tested by an accredited testing institute and complies with the requirements for Medical Device Directive Class I products (MDD 93/42/EEC). Sabina™ 200 complies with the requirements for IEC 60601-1, IEC 60601-1-2, EN ISO 10535, UL-60601-1 and CAN/CSA C22.2 No.601.1.

Design and Quality by Liko in Sweden.

Liko is quality-certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485. Liko is also certified in accordance with environmental standard ISO 14001.

Hill-Rom and Design and Hill-Rom word mark are trademarks or registered trademarks of Hill-Rom Services Inc. Liko and Design and Liko word mark are trademarks or registered trademarks of Liko R&D AB. Sabina, SupportVest and SafetyVest are trademarks or registered trademarks of Liko R&D AB.



Enhancing Outcomes for Patients and Their Caregivers.™

- France \_\_\_\_\_ +33 (0)2 97 50 92 12
- United Kingdom \_\_\_\_\_ +44 (0)1530 411000
- Deutschland \_\_\_\_\_ +49 (0)211 16450-0
- Nederland \_\_\_\_\_ +31 (0)347 32 35 32
- Italia \_\_\_\_\_ + 39 (0)2 950541
- Suisse/Schweiz \_\_\_\_\_ +41 (0)21 706 21 30
- (deutschsprachig) \_\_\_\_\_ +41 (0)21 706 21 38

- Österreich \_\_\_\_\_ +43 (0)2243 28550
- Ireland \_\_\_\_\_ +353 (0)1 413 6005
- Iberia \_\_\_\_\_ +34 (0)93 6856000
- Nordic Countries \_\_\_\_\_ +46 (0)920 47 47 11
- Export \_\_\_\_\_ Contact your local Hill-Rom distributor or contact your Area Manager via website or call \_\_\_\_\_ +1 812 934 8173

[www.hill-rom.com](http://www.hill-rom.com)  
[www.liko.com](http://www.liko.com)

Autorizované zastoupení pro ČR:



**DARTIN, spol. s r.o.**  
Suchbolská 688, 25262 Horoměřice  
Tel./fax: 241 470 361, 261 216 076  
[www.dartin.cz](http://www.dartin.cz)  
e-mail: [info@dartin.cz](mailto:info@dartin.cz)

Autorizované zastúpenie pre SR:



**DARTIN Slovensko, spol. s r.o.**  
L. Podjavorinskej 25, 08005 Prešov  
Tel./fax: 051/7734403, 7723634  
[www.dartin.sk](http://www.dartin.sk)  
e-mail: [info@dartin.sk](mailto:info@dartin.sk)