


High Flow Therapy

# Take Comfort

High Flow Therapy  
through a Single Prong



## SOLO™ Cannula




**A whole new way to assure an open system  
while delivering High Flow Therapy**

- Assures an open system
- Eliminates concerns about over occlusion
- Simplifies NG tube placement
- Assures less than 50% occlusion in even the smallest baby
- May facilitate high flow therapy in patients with anatomical defects

## High Flow Therapy through a Single Prong Cannula

SOLO™ works with the Vapotherm Precision Flow to provide the same gentle and effective ventilation support of dual prong cannulas. Only Precision Flow can deliver High Flow Therapy through a single prong.

- Avoid masks and associated skin breakdown
- Avoid intubation and associated costs and adverse events
- Easily access, feed and hold the patient



“The single prong cannula has worked very well for us particularly to support neonates under 1000 grams. With the single prong cannula we do not worry about occluding too much of the nares.”

– **Jorge Rojas, MD**  
Staff Neonatologist

### Technical Specifications

|               |                |
|---------------|----------------|
| Gas flow rate | 0 – 8 lpm      |
| Tip diameter  | 0.073" (1.9mm) |

### Ordering Information

|  |                 |
|--|-----------------|
| SOLO Single Prong Infant / Neonate / Premature Nasal Cannula | <b>SOLO1300</b> |
|--|-----------------|

Autorizované zastoupení pro ČR:

**DARTiN**

DARTIN, spol. s r.o.  
Suchbolská 688, 252 62 Horoměřice  
Tel./fax: 241 470 361, 261 216 076  
www.dartin.cz  
e-mail: info@dartin.cz

Autorizované zastúpenie pre SR:

**DARTiN**

DARTIN Slovensko, s.r.o.  
L. Podjavorinskej 25, 080 05 Prešov  
Tel./fax: 051/7734403, 7723634  
www.dartin.sk  
e-mail: info@dartin.sk



**VAPOTHERM®**

Vapotherm, Inc. • 22 Industrial Drive • Exeter, NH 03833 • 866.827.6843 • www.vtherm.com

3100185 Rev A 7/13